

CENTRAL VENOUS CATHETER

(Available at Present)

P - Proximal, D - Distal, M-Medial

SINGLE LUMEN CENTRAL VENOUS CATHETER

08.01.2009

OD (FR) / Ga (I.D.)	Length (mm)	GUIDEWIRE		Needle	Dilator OD (FR) / Length (mm)	Syringe	Qty/Box
		OD (inch/mm)	Length (mm)				
5FR / 16G (D)	160mm	0.032 / 0.81	45 / 450 & 50/500	18g	6FR / 82mm	5cc	10
5FR / 16G (D)	160mm	0.035 / 0.89	45 / 450 & 50/500	18g	6FR / 82mm	5cc	10
5FR / 16G (D)	200mm	0.032 / 0.81	45 / 450 & 50/500	18g	6FR / 82mm	5cc	10
5FR / 16G (D)	200mm	0.035 / 0.89	45 / 450 & 50/500	18g	6FR / 82mm	5cc	10
5FR / 16G (D)	300mm	0.032 / 0.81	70/700mm	18g	6FR / 82mm	5cc	10
5FR / 16G (D)	300mm	0.035 / 0.89	70/700mm	18g	6FR / 82mm	5cc	10

DOUBLE LUMEN CENTRAL VENOUS CATHETER

OD (FR) / Ga (I.D.)	Length (mm)	GUIDEWIRE		Needle	Dilator OD (FR) / Length (mm)	Syringe	Qty/Box
		OD (inch/mm)	Length (mm)				
7FR / 16G(P) - 16G(D)	160mm	0.032 / 0.81	45 / 450 & 50/500	18g	8FR / 100mm	5cc	10
7FR / 16G(P) - 16G(D)	160mm	0.035 / 0.89	45 / 450 & 50/500	18g	8FR / 100mm	5cc	10
7FR / 16G(P) - 16G(D)	200mm	0.032 / 0.81	45 / 450 & 50/500	18g	8FR / 100mm	5cc	10
7FR / 16G(P) - 16G(D)	200mm	0.035 / 0.89	45 / 450 & 50/500	18g	8FR / 100mm	5cc	10
7FR / 16G(P) - 16G(D)	300mm	0.032 / 0.81	70/700mm	18g	8FR / 100mm	5cc	10
7FR / 16G(P) - 16G(D)	300mm	0.035 / 0.89	70/700mm	18g	8FR / 100mm	5cc	10

TRIPLE LUMEN CENTRAL VENOUS CATHETER

OD (FR) / Ga (I.D.)	Length (mm)	GUIDEWIRE		Needle	Dilator OD (FR) / Length (mm)	Syringe	Qty/Box
		OD (inch/mm)	Length (mm)				
7FR / 18G(P)/16G(D)/18G(M)	160mm	0.032 / 0.81	45 / 450 & 50/500	18g	8FR / 100mm	5cc	10
7FR / 18G(P)/16G(D)/18G(M)	160mm	0.035 / 0.89	45 / 450 & 50/500	18g	8FR / 100mm	5cc	10
7FR / 18G(P)/16G(D)/18G(M)	200mm	0.032 / 0.81	45 / 450 & 50/500	18g	8FR / 100mm	5cc	10
7FR / 18G(P)/16G(D)/18G(M)	200mm	0.035 / 0.89	45 / 450 & 50/500	18g	8FR / 100mm	5cc	10
7FR / 18G(P)/16G(D)/18G(M)	300mm	0.032 / 0.81	70/700mm	18g	8FR / 100mm	5cc	10
7FR / 18G(P)/16G(D)/18G(M)	300mm	0.035 / 0.89	70/700mm	18g	8FR / 100mm	5cc	10

CENTRAL VENOUS CATHETER
(Available by March)

P - Proximal, D - Distal, M-Medial

SINGLE LUMEN CENTRAL VENOUS CATHETER

OD (FR) / Ga (I.D.)	Length (mm)	GUIDEWIRE		Needle	Dilator OD (FR) / Length (mm)	Syringe	Qty/Box
		OD (inch/mm)	Length (mm)				
6FR / 14G (D)	160mm	0.032 / 0.81	45 / 450	18g	7FR / 82mm	5cc	10
6FR / 14G (D)	160mm	0.035 / 0.89	45 / 450	18g	7FR / 82mm	5cc	10
6FR / 14G (D)	200mm	0.032 / 0.81	45 / 450	18g	7FR / 82mm	5cc	10
6FR / 14G (D)	200mm	0.035 / 0.89	45 / 450	18g	7FR / 82mm	5cc	10
6FR / 14G (D)	300mm	0.032 / 0.81	70/700mm	18g	7FR / 82mm	5cc	10
6FR / 14G (D)	300mm	0.035 / 0.89	70/700mm	18g	7FR / 82mm	5cc	10

DOUBLE LUMEN CENTRAL VENOUS CATHETER

OD (FR) / Ga (I.D.)	Length (mm)	GUIDEWIRE		Needle	Dilator OD (FR) / Length (mm)	Syringe	Qty/Box
		OD (inch/mm)	Length (mm)				
7FR / 14G(P) - 18G(D)	160mm	0.025 / 0.64	45 / 450	18g	8FR / 100mm	5cc	10
7FR / 14G(P) - 18G(D)	200mm	0.025 / 0.64	45 / 450	18g	8FR / 100mm	5cc	10
7FR / 14G(P) - 18G(D)	300mm	0.025 / 0.64	70/700mm	18g	8FR / 100mm	5cc	10

TRIPLE LUMEN CENTRAL VENOUS CATHETER

OD (FR) / Ga (I.D.)	Length (mm)	GUIDEWIRE		Needle	Dilator OD (FR) / Length (mm)	Syringe	Qty/Box
		OD (inch/mm)	Length (mm)				
5.5FR / 22G(P)/20G(D)/22G(M)	80mm	0.018/0.46	45 / 450	21g	6FR / 52mm	5cc	10
5.5FR / 22G(P)/20G(D)/22G(M)	160mm	0.018/0.46	45 / 450	21g	6FR / 52mm	5cc	10

REMARKS : All CVC are with Y HUB with NRV & BD Syringe