

## Art-Cath

| Ref Code | Fr | Size | Catheter Length (mm) | Flowrate (ml/min) | Guidewire Length (mm) | Introducer Needle | Packaging Duplex / Carton |
|----------|----|------|----------------------|-------------------|-----------------------|-------------------|---------------------------|
| 13528    | 3  | 20G  | 40                   | 38                | 200                   | 20Gx38mm          | 20/200                    |
| 13371    | 3  | 20G  | 60                   | 27                | 200                   | 20Gx38mm          | 20/200                    |
| 13529    | 3  | 20G  | 80                   | 33                | 300                   | 20Gx38mm          | 20/200                    |
| 13372    | 4  | 18G  | 100                  | 24                | 300                   | 19Gx52mm          | 20/200                    |
| 13534    | 4  | 18G  | 180                  | 95                | 460                   | 18Gx70mm          | 20/200                    |
| 13410    | 5  | 16G  | 150                  | 29                | 390                   | 18Gx70mm          | 20/200                    |