

Autofusion Nuovo Set		
Product Code	Product Description	Units / Duplex Box
14263	Autofusion Nuovo IV Infusion Set, Vented Spike, Extruded drip chamber 63mm with autostop filter, DEHP free tube 180cm, Orange Roller body, Rotating luer lock with PE Cap -Soft Blister pack	30
14279	AUTOFUSION NUOVO IV fluid Infusion set with Vented spike, 180cm DEHP free tube, with Single Drip Chamber with Air Stop filter with long roller clamp, without Rotating Luer Lock with Hydrophobic Filter cap- Paper Pack	30
14293	Autofusion Nuovo IV Infusion Set, Vented Spike, Extruded drip chamber 63mm with autostop filter, DEHP free tube 180cm, Orange Roller body, Latex free Y-site, Rotating luer lock with PE Cap -Soft Blister pack	30
14294	Autofusion Nuovo IV Infusion Set, Vented Spike, Extruded drip chamber 63mm with autostop filter, DEHP free tube 180cm, Orange Roller body, Needle free Y-site, Rotating luer lock with PE Cap -Soft Blister pack	30
14296	Autofusion Nuovo IV Infusion Set, Vented Spike, Extruded drip chamber 63mm with Autostop filter, DEHP free tube 180cm, Orange Roller body, Rotating luer lock with PE Cap, Latex free Y-site -Soft Blister pack	30