

Autofusion Set		
Product Code	Product Description	Units / Duplex Box
14260	Auto stop IV Infusion Set, with transparent Vented spike, PVC free double Drip Chamber with Autostop filter, DEHP free tube 150cm, Long White Roller Clamp, Rotating Luer Lock with Purge Filter Cap- Blister Pack	30
14261	Autofusion IV Infusion Set, Vented spike, PVC free double Drip Chamber with Autostop filter, DEHP free tube 180cm, Latex free Y-site, Long Orange Roller Clamp, Rotating Luer Lock with PE cap with filter- Blister Pack	30
14262	Autofusion IV Infusion Set, Vented Spike, PVC Free double drip chamber with autostop filter, DEHP free tube 180cm, Long Orange Roller Clamp, BACK CHECK VALVE with PE Cap with filter-Soft Blister pack	30
14263	Autofusion Nuovo IV Infusion Set, Vented Spike, Extruded drip chamber 63mm with autostop filter, DEHP free tube 180cm, Orange Roller body, Rotating luer lock with PE Cap -Soft Blister pack	30
14264	Autofusion IV Infusion Set, Vented Spike, PVC Free double drip chamber with autostop filter, Needle free Y-site, DEHP free tube 210cm, Long Orange Roller Clamp, Rotating luer lock with Purge filter- Blister pack	30
14268	Autofusion IV Infusion Set, Vented Spike, PVC Free double drip chamber with autostop filter, Three way stop cock, DEHP free tube 220cm, Long Orange Roller Clamp, Rotating luer lock with Purge filter- Blister pack	30
14269	Autofusion IV Infusion Set, Vented Spike, PVC Free double drip chamber with autostop filter, Needle free Y-site, DEHP free tube 200cm, Long Orange Roller Clamp, Rotating luer lock with Purge filter- Soft Blister pack	30