

Oncofusion Set		
Product Code	Product Description	Units / Duplex Box
14160	IV Infusion Set, 67 mm drip chamber,Vented, 0.2µm Filter, Needle free y-site, White roller clamp, Rotating Luer Lock, PE cap with filter, PVC free Tube-150cm-Blister Pack	10
14161	IV Infusion Set, 67 mm drip chamber,Vented, 0.2µm Filter, White roller clamp, Rotating Luer Lock, PE cap with filter, DEHP free Tube-150cm-Blister Pack	10
14162	IV Infusion Set, 67 mm drip chamber,Vented, 0.2µm Filter, White roller clamp, Rotating Luer Lock, PE cap with filter, PVC free Tube-150cm-Blister Pack	10
14163	IV Infusion Set, 67 mm drip chamber, Vented, 0.2µm Filter, Latex free y-site, White roller clamp, Rotating Luer Lock, PE cap with filter, PVC free Tube-150cm-Blister Pack	10
14164	IV Infusion Set , with transparent Vented spike, PVC free double Drip Chamber with Autostop filter, 0.2µm inline Filter, White roller clamp, fixed Luer Lock, PE cap with filter, PVC DEHP free Tube-120+30cm-Blister Pack	10