

Polycath IV Cannula						
Ref. Code	Gauge	Catheter Length (mm)	Catheter OD (mm)	Flow Rate (ml/min)	Catheter Material	Packaging
10410	16 	45	1.70	200	FEP	100 each/Duplex & 10 Duplex/Case
10411	17 	45	1.50	142	FEP	100 each/Duplex & 10 Duplex/Case
10412	18 	45	1.30	90	FEP	100 each/Duplex & 10 Duplex/Case
10413	18 	32	1.30	90	FEP	100 each/Duplex & 10 Duplex/Case
10414	20 	32	1.10	61	FEP	100 each/Duplex & 10 Duplex/Case
10415	22 	25	0.90	36	FEP	100 each/Duplex & 10 Duplex/Case
10420	16 	45	1.70	200	PUR	100 each/Duplex & 10 Duplex/Case
10421	17 	45	1.50	142	PUR	100 each/Duplex & 10 Duplex/Case
10422	18 	45	1.30	90	PUR	100 each/Duplex & 10 Duplex/Case
10423	18 	32	1.30	90	PUR	100 each/Duplex & 10 Duplex/Case
10424	20 	32	1.10	61	PUR	100 each/Duplex & 10 Duplex/Case
10425	22 	25	0.90	36	PUR	100 each/Duplex & 10 Duplex/Case

Packaging Dimension	
Unit Pack Size	114mm X 27mm
Duplex Box Size	170mm X 117mm X 150mm
Case Size	605mm X 320mm X 188mm