

Poly Lead-OTW		
Product Code	Product Description	Units/Box
13535	Mini midline over the Wire FR-4, Size-18G, catheter-40mm, Guidewire length-250mm, Needle-20Gx54mm	100
13536	13536 Mini midline over the Wire FR-4, Size-18G, catheter-60mm, Guidewire length-250mm, Needle-20Gx54mm	100
13537	Mini midline over the Wire FR-4, Size-18G, catheter-80mm, Guidewire length-250mm, Needle-20Gx54mm	100
13538	Mini midline over the Wire FR-4, Size-18G, catheter-120mm, Guidewire length-400mm, Needle-20Gx54mm	100
13539	Mini midline over the Wire FR-3, Size-20G, catheter-40mm, Guidewire length-250mm, Needle-20Gx54mm	100
13580	Mini midline over the Wire FR-3, Size-20G, catheter-60mm, Guidewire length-250mm, Needle-20Gx54mm	100
13581	Mini midline over the Wire FR-3, Size-20G, catheter-80mm, Guidewire length-250mm, Needle-20Gx54mm	100
13582	Mini midline over the Wire FR-3, Size-20G, catheter-120mm, Guidewire length-400mm, Needle-20Gx54mm	100
13583	Mini midline over the Wire FR-2, Size-22G, catheter-40mm, Guidewire length-250mm, Needle-20Gx42mm	100
13584	Mini midline over the Wire FR-2, Size-22G, catheter-60mm, Guidewire length-250mm, Needle-20Gx42mm	100
13585	Mini midline over the Wire FR-2, Size-22G, catheter-80mm, Guidewire length-250mm, Needle-20Gx42mm	100
13586	Mini midline over the Wire FR-2, Size-22G, catheter-120mm, Guidewire length-400mm, Needle-20Gx42mm	100