

Polypen IV Cannula						
Ref. Code	Gauge	Catheter Length (mm)	Catheter OD (mm)	Flow Rate (ml/min)	Catheter Material	Packaging
10290	14 	45	2.10	305	FEP	100 each/Duplex & 10 Duplex/Case
10291	16 	45	1.70	200	FEP	100 each/Duplex & 10 Duplex/Case
10292	17 	45	1.50	142	FEP	100 each/Duplex & 10 Duplex/Case
10293	18 	45	1.30	90	FEP	100 each/Duplex & 10 Duplex/Case
10294	18 	32	1.30	90	FEP	100 each/Duplex & 10 Duplex/Case
10295	20 	32	1.10	61	FEP	100 each/Duplex & 10 Duplex/Case
10296	22 	25	0.90	36	FEP	100 each/Duplex & 10 Duplex/Case
10297	24 	19	0.70	23	FEP	100 each/Duplex & 10 Duplex/Case
10298	26 	19	0.60	17	FEP	100 each/Duplex & 10 Duplex/Case
10640	14 	45	2.10	305	PUR	100 each/Duplex & 10 Duplex/Case
10641	16 	45	1.70	200	PUR	100 each/Duplex & 10 Duplex/Case
10642	17 	45	1.50	142	PUR	100 each/Duplex & 10 Duplex/Case
10643	18 	45	1.30	90	PUR	100 each/Duplex & 10 Duplex/Case
10644	18 	32	1.30	90	PUR	100 each/Duplex & 10 Duplex/Case
10645	20 	32	1.10	61	PUR	100 each/Duplex & 10 Duplex/Case
10646	22 	25	0.90	36	PUR	100 each/Duplex & 10 Duplex/Case
10647	24 	19	0.70	23	PUR	100 each/Duplex & 10 Duplex/Case
10648	26 	19	0.60	17	PUR	100 each/Duplex & 10 Duplex/Case

Packaging Dimension	
Unit Pack Size	114mm X 27mm
Duplex Box Size	170mm X 117mm X 150mm
Case Size	605mm X 320mm X 188mm