

Polytrol - Micro Set		
Product Code	Product Description	Units / Duplex Box
14280	IV Infusion Set, Vented Spike with Micropin needle, 57mm Drip Chamber, Flow regulator, White click clamp, Thumb Y-injection port, Luer Lock, DEHP Free Tube 150cm, Blister Pack	25
14281	IV Infusion Set, Non-Vented Spike with Micropin needle, 57mm Drip Chamber, Flow regulator, Blue click clamp, Latex free Y injection port, Luer Lock, DEHP Free Tube 150cm, Blister Pack	25
14282	IV Infusion Set, Non-Vented Spike, 57mm Drip Chamber, Flow regulator, Blue click clamp, Luer Lock, DEHP Free Tube 150cm, Blister Pack	25