

SC-Compliance Chart

Inflation Pressure		Balloon Diameter (mm)									
kPa	(atm)	1.50 mm	2.00 mm	2.25 mm	2.50 mm	2.75 mm	3.00 mm	3.25 mm	3.50 mm	3.75 mm	4.00 mm
405	4	1.42	1.91	2.18	2.41	2.66	2.91	3.15	3.42	3.63	3.88
507	5	1.46	1.95	2.22	2.45	2.71	2.96	3.21	3.46	3.68	3.94
NP*608	6	1.51	1.99	2.26	2.49	2.77	3.02	3.26	3.53	3.73	4.02
709	7	1.55	2.03	2.3	2.53	2.81	3.07	3.31	3.58	3.78	4.07
811	8	1.59	2.07	2.34	2.57	2.85	3.12	3.37	3.63	3.82	4.12
912	9	1.63	2.11	2.38	2.61	2.89	3.17	3.41	3.68	3.87	4.18
1013	10	1.67	2.15	2.41	2.66	2.93	3.21	3.46	3.72	3.91	4.24
1115	11	1.71	2.19	2.45	2.7	2.97	3.26	3.51	3.77	3.96	4.3
1216	12	1.75	2.23	2.49	2.75	3.01	3.31	3.56	3.81	4.01	4.35
1317	13	1.79	2.27	2.53	2.79	3.05	3.36	3.6	3.86	4.06	4.41
RBP**1419	14	1.83	2.31	2.57	2.83	3.08	3.41	3.64	3.91	4.12	4.46
1520	15	1.88	2.35	2.61	2.88	3.13	3.45	3.69	3.95	4.17	4.51
1621	16	1.92	2.39	2.66	2.92	3.17	3.5	3.74	3.99	4.22	4.56
1723	17	1.96	2.43	2.7	2.96	3.21	3.54	3.78	4.04	4.26	4.61

