

Stent Compliance Chart

Inflation Pressure	Stent Diameter Ø						
(atm/bar)	2.00	2.25	2.50	2.75	3.00	3.50	4.00
6.00	1.80	2.10	2.35	2.62	2.83	3.42	3.80
*8.00	2.00	2.25	2.48	2.74	3.04	3.52	4.02
10.00	2.15	2.45	2.60	2.81	3.12	3.62	4.10
12.00	2.20	2.53	2.70	2.95	3.24	3.85	4.25
14.00	2.25	2.60	2.80	3.05	3.38	4.00	4.35
*16.00	2.32	2.70	2.90	3.15	3.50	4.08	4.49
18.00	2.37	2.75	3.00	3.20	3.60	4.13	4.69
20.00	2.42	2.77	3.05	3.27	3.65	4.22	4.72
22.00	2.45	2.80	3.08	3.32	3.68	4.30	4.75