

NC-Compliance Chart

Inflation Pressure		Balloon Diameter (mm)									
kPa	(atm)	1.50 mm	2.00 mm	2.25 mm	2.50 mm	2.75 mm	3.00 mm	3.25 mm	3.50 mm	3.75 mm	4.00 mm
811	8	1.44	1.96	2.15	2.41	2.71	2.93	3.18	3.42	3.68	3.92
912	9	1.47	1.99	2.19	2.47	2.74	2.96	3.21	3.45	3.72	3.96
NP*1013	10	1.5	2.03	2.22	2.52	2.78	3	3.25	3.48	3.76	4.01
1115	11	1.52	2.05	2.25	2.55	2.81	3.03	3.28	3.52	3.8	4.05
1216	12	1.54	2.07	2.27	2.57	2.83	3.06	3.32	3.55	3.83	4.08
1317	13	1.56	2.09	2.29	2.6	2.86	3.09	3.35	3.58	3.86	4.12
1419	14	1.58	2.11	2.31	2.63	2.88	3.12	3.38	3.62	3.9	4.16
1520	15	1.6	2.13	2.33	2.65	2.9	3.16	3.41	3.66	3.93	4.2
1621	16	1.62	2.15	2.35	2.68	2.92	3.19	3.45	3.69	3.97	4.23
1723	17	1.64	2.17	2.37	2.71	2.94	3.22	3.48	3.72	4	4.31
RBP**1824	18	1.66	2.19	2.39	2.74	2.96	3.25	3.51	3.76	4.04	4.35
1925	19	1.68	2.21	2.41	2.78	2.99	3.28	3.55	3.8	4.08	4.38
2027	20	1.7	2.23	2.43	2.82	3.01	3.32	3.59	3.83	4.11	4.41
2128	21	1.72	2.25	2.45	2.85	3.04	3.36	3.63	3.86	4.15	4.45
2229	22	1.74	2.27	2.47	2.87	3.06	3.4	3.66	3.9	4.17	4.49

