

NEOFOLD PRO

Manifold

The Manifold features a clear body for visibility, easy rotating handles with balanced pressure and rotation, a standard connection adaptor compatible with various devices, and 60°C & 180°C rotating handles for enhanced rotatability. It is with pressure rating up to 250 and 500 PSI.

Available in:
2 Port Right ON, 3 Port Right ON, 2 Port Right OFF, 3 Port Right OFF



| Neofold Pro Ordering Information | | | | |
|----------------------------------|--------------|----------|-----------|------------|
| Product | Product Code | Pressure | Side Port | Valve Type |
| Manifold | 13755 | 250psi | 2 | OFF |
| | 13753 | 250psi | 2 | ON |
| | 13756 | 250psi | 3 | OFF |
| | 13754 | 250psi | 3 | ON |
| | 13757 | 500psi | 2 | OFF |
| | 13759 | 500psi | 2 | ON |
| | 13758 | 500psi | 3 | OFF |
| | 13760 | 500psi | 3 | ON |