

Technical Information Excedo SC	
Balloon Diameter (mm)	1.25, 1.50, 2.00, 2.25, 2.75, 3.00, 3.25, 3.50, 3.75, 4.00
Balloon Length (mm)	6, 8, 10, 12, 15, 18, 20, 21, 25, 30
Crossing Profile	0.027" (1.50mm)
Tip Entry Profile	0.019"
Nominal Pressure	06 atm
Rated Burst Pressure	14 atm
Balloon Material	Pebax
Marker Bands	Two Gold Markers
Folding	3-fold
Coating	Hydrophilic
Guiding Catheter Compatibility	6F
Guide wire Diameter	0.014"
Proximal shaft Diameter OD	1.81FR
Distal Shaft Diameter OD	2.26FR
Usable catheter length (cm)	145
Tapered tip	Yes

