

Technical Information Excedo NC	
Balloon Diameter (mm)	1.50, 2.00, 2.25, 2.75, 3.00, 3.25, 3.50, 3.75, 4.00
Balloon Length (mm)	6, 8, 9, 10, 12,15,18,19, 20, 21, 25, 27
Crossing Profile	0.027" (1.50mm)
Tip Entry Profile	0.019"
Nominal Pressure	10 atm
Rated Burst Pressure	18 atm
Balloon Material	Nylon
Marker Bands	Two Gold Markers
Folding	3-fold
Coating	Hydrophilic
Guiding Catheter Compatibility	6F
Guide wire Diameter	0.014"
Proximal shaft Diameter OD	1.81FR
Distal Shaft Diameter OD	2.26FR
Usable catheter length (cm)	145
Tapered tip	Yes