

TERAVOX


Inflation Device

Features:

- Confident-grip torque knob with a ridged design for “no slip” torque control and accurate pressure adjustment
- Ergonomic pistol grip with centered trigger that fits comfortably in hand
- Luminescent dial appearance that provides proper readability in dim light conditions
- Superb visibility with black & white dial pressure gauge provides reading support
- Angled pressure gauge provides a suitable reading angle for operating procedures

Available in:

30 atm pressure gauge
with 20 ml clear syringe barrel



Teravox Ordering Information		
Product	Product Code	Variant
Inflation Device	90284	Plunger type
	88733	Gun type