

Central Venous Catheters Sizes

PRODUCT DESCRIPTION		INDWELLING CATHETER		GUIDE WIRE	
Variant	Novocent	Specification (Sizes)	Length (mm)	Diameter (inch/mm)	Length (cm/mm)
Single Lumen	15001	3FR / 20G (D)	80	0.018 / 0.46	50/500
Single Lumen	15002	3FR / 20G (D)	130	0.018 / 0.46	50/500
Single Lumen	15003	3FR / 20G (D)	150	0.018 / 0.46	50/500
Single Lumen	15004	3FR / 20G (D)	160	0.018 / 0.46	50/500
Single Lumen	15005	3FR / 20G (D)	200	0.018 / 0.46	50/500
Single Lumen	15006	4FR / 18G (D)	130	0.018 / 0.46	50/500
Single Lumen	15007	4FR / 18G (D)	150	0.018 / 0.46	50/500
Single Lumen	15008	4FR / 18G (D)	160	0.018 / 0.46	50/500
Single Lumen	15009	4FR / 18G (D)	200	0.018 / 0.46	50/500
Single Lumen	15010	5FR / 16G (D)	160	0.035 / 0.89	50/500
Single Lumen	15011	5FR / 16G (D)	200	0.035 / 0.89	50/500
Single Lumen	15012	5FR / 16G (D)	300	0.035 / 0.89	70/700
Single Lumen	15013	6FR / 14G (D)	160	0.035 / 0.89	50/500
Single Lumen	15014	6FR / 14G (D)	200	0.035 / 0.89	50/500
Single Lumen	15015	6FR / 14G (D)	300	0.035 / 0.89	70/700
Single Lumen	15016	7FR / 12G (D)	160	0.035 / 0.89	50/500
Single Lumen	15017	7FR / 12G (D)	200	0.035 / 0.89	50/500
Single Lumen	15019	8FR / 12G (D)	200	0.035 / 0.89	50/500
Single Lumen	15020	9FR / 11G (D)	200	0.035 / 0.89	50/500

Packaging Dimension	
Unit Pack Size	235mm X 170mm
Duplex Box Size	245mm X240mm X 180mm
Case Size	950mm X 510mm X 265mm

Central Venous Catheters Sizes

PRODUCT DESCRIPTION		INDWELLING CATHETER		GUIDE WIRE	
Variant	Novocent	Specification (Sizes)	Length (mm)	Diameter (inch/mm)	Length (cm/mm)
Double Lumen	15068	4FR / 22G (P) - 22 G (D)	80	0.018 / 0.46	50/500
Double Lumen	15050	4FR / 22G (P) - 22 G (D)	130	0.018 / 0.46	50/500
Double Lumen	15051	4FR / 22G (P) - 22 G (D)	160	0.018 / 0.46	50/500
Double Lumen	15052	4FR / 22G (P) - 22 G (D)	200	0.018 / 0.46	50/500
Double Lumen	15069	5FR / 20G (P) - 18 G (D)	80	0.018 / 0.46	50/500
Double Lumen	15053	5FR / 20G (P) - 18 G (D)	130	0.018 / 0.46	50/500
Double Lumen	15064	5FR / 20G (P) - 18 G (D)	160	0.018 / 0.46	50/500
Double Lumen	15054	5FR / 20G (P) - 18 G (D)	200	0.018 / 0.46	50/500
Double Lumen	15055	7FR / 16G (P) - 16 G (D)	160	0.035 / 0.89	50/500
Double Lumen	15056	7FR / 16G (P) - 16 G (D)	200	0.035 / 0.89	50/500
Double Lumen	15057	7FR / 16G (P) - 16 G (D)	300	0.035 / 0.89	70/700
Double Lumen	15058	7FR / 14G (P) - 18 G (D)	160	0.025 / 0.64	50/500
Double Lumen	15059	7FR / 14G (P) - 18 G (D)	200	0.025 / 0.64	50/500
Double Lumen	15060	7FR / 14G (P) - 18 G (D)	300	0.025 / 0.64	70/700
Double Lumen	15065	7FR / 18G (P) - 14 G (D)	160	0.035 / 0.89	50/500
Double Lumen	15066	7FR / 18G (P) - 14 G (D)	200	0.035 / 0.89	50/500
Double Lumen	15061	8FR / 14G (P) - 14 G (D)	160	0.035 / 0.89	50/500
Double Lumen	15063	8FR / 14G (P) - 14 G (D)	300	0.035 / 0.89	70/700

Packaging Dimension	
Unit Pack Size	235mm X 170mm
Duplex Box Size	245mm X240mm X 180mm
Case Size	950mm X 510mm X 265mm

Central Venous Catheters Sizes

PRODUCT DESCRIPTION		INDWELLING CATHETER		GUIDE WIRE	
Variant	Novocent	Specification (Sizes)	Length (mm)	Diameter (inch/mm)	Length (cm/mm)
Triple Lumen	15106	4.5FR / 23G (P) 20G (D)-23G(M)	60	0.018 / 0.46	50/500
Triple Lumen	15115	4.5FR / 23G (P) 20G (D)-23G(M)	80	0.018 / 0.46	50/500
Triple Lumen	15116	4.5FR / 23G (P) 20G (D)-23G(M)	130	0.018 / 0.46	50/500
Triple Lumen	15108	5FR / 22G (P) 20G (D)-22G(M)	80	0.018 / 0.46	50/500
Triple Lumen	15109	5FR / 22G (P) 20G (D)-22G(M)	130	0.018 / 0.46	50/500
Triple Lumen	15110	5FR / 22G (P) 20G (D)-22G(M)	150	0.018 / 0.46	50/500
Triple Lumen	15111	5FR / 22G (P) 20G (D)-22G(M)	160	0.018 / 0.46	50/500
Triple Lumen	15112	5FR / 22G (P) 20G (D)-22G(M)	200	0.018 / 0.46	50/500
Triple Lumen	15100	5.5FR / 22G (P) 20G (D)-22G(M)	80	0.018 / 0.46	50/500
Triple Lumen	15101	5.5FR / 22G (P) 20G (D)-22G(M)	160	0.018 / 0.46	50/500
Triple Lumen	15102	5.5FR / 22G (P) 20G (D)-22G(M)	200	0.018 / 0.46	50/500
Triple Lumen	15113	7FR / 18G (P) 16G (D)-18G(M)	150	0.035 / 0.89	50/500
Triple Lumen	15103	7FR / 18G (P) 16G (D)-18G(M)	160	0.035 / 0.89	50/500
Triple Lumen	15104	7FR / 18G (P) 16G (D)-18G(M)	200	0.035 / 0.89	50/500
Triple Lumen	15105	7FR / 18G (P) 16G (D)-18G(M)	300	0.035 / 0.89	70/700
Four Lumen	-	8FR / 18G (P) 16G (D)-18G(M1-M2)	130	0.035 / 0.89	50/500
Four Lumen	-	8FR / 18G (P) 16G (D)-18G(M1-M2)	160	0.035 / 0.89	50/500
Four Lumen	-	8FR / 18G (P) 16G (D)-18G(M1-M2)	200	0.035 / 0.89	50/500

Packaging Dimension	
Unit Pack Size	235mm X 170mm
Duplex Box Size	245mm X240mm X 180mm
Case Size	950mm X 510mm X 265mm

Central Venous Catheters Sizes

Sr No.	Generic Name of Device	Brand Name of Device	Product Code
1	Central Venous Catheter Single Lumen	Novocent	15001-15025, 15031-15049, 15350
2	Central Venous Catheter Single Lumen	Novocent protekt	17000-17023, 18050-18064
3	Central Venous Catheter Double Lumen	Novocent	15050-15069, 15071, 15085-15099, 15361-15363
4	Central Venous Catheter Double Lumen	Novocent protekt	17030-17050, 18070-18082
5	Central Venous Catheter Triple Lumen	Novocent	15100-15117, 15150-15163
6	Central Venous Catheter Triple Lumen	Novocent protekt	17055-17073, 18085-18093
7	Central Venous Catheter Quintuple Lumen	Novocent	17271-17272
8	Central Venous Catheter Quintuple Lumen	Novocent protekt	17260-17270, 18100-18102, 18110-18112

Packaging Dimension	
Unit Pack Size	235mm X 170mm
Duplex Box Size	245mm X240mm X 180mm
Case Size	950mm X 510mm X 265mm